



































































Replacement overview

previous hearing aid family	new hearing aid family	new hearing aid family overview (models available in each Avio series)		controls (standard/option)	gain [dB]	max out-put [dB]	fitting range
SHAPE	avio ⁵	 CIC 10; CIC 10 HPG	A510 / A510 HPG A310 / A310 HPG A110 / A110 HPG	on/off battery door	41 / 49	118 / 123	
			A510 PB / A510 PB HPG A310 PB / A310 PB HPG A110 PB / A110 PB HPG	on/off battery door; push button	41 / 50	119 / 124	
		 ITC 30; ITC 30 HPG	A530 / A530 HPG A330 / A330 HPG A130 / A130 HPG	on/off battery door; push button; VC wheel	52 / 58	121 / 128	
			A530-D / A530-D HPG A330-D / A330-D HPG A130-D / A130-D HPG	on/off battery door; push button; VC wheel	52 / 58	121 / 128	
		 ITE HS 40; ITE HF 40 HPG	A540 / A540 HPG A340 / A340 HPG A140 / A140 HPG	on/off battery door; push button; VC wheel	61 / 67	129 / 134	
			A540-D / A540-D HPG A340-D / A340-D HPG A140-D / A140-D HPG	on/off battery door; push button; VC wheel	61 / 67	129 / 134	
SLIM	avio ³	 ITE FS 50; ITE FS 50 HPG	A550 / A550 HPG A350 / A350 HPG A150 / A150 HPG	on/off battery door; VC wheel; push button;	61 / 67	129 / 134	
			A550-D / A550-D HPG A350-D / A350-D HPG A150-D / A150-D HPG	on/off battery door; push button; VC wheel	61 / 67	129 / 134	
		 RIE	A560-D / A360-D / A160-D	on/off battery door; push button; VC toggle	64	127	
SNAP	avio ³	 RIE open	A560-D open / A360-D open / A160-D open	on/off battery door; push button; VC toggle	61	125	
		 RIE HPG	A560-D HPG / A360-D HPG / A160-D HPG	on/off battery door; push button; VC toggle	68	130	
IQ	avio ¹	 mini BTE	A565-D / A365-D / A165-D	on/off battery door; push button	62	132	
		 mini BTE open	A565-D open / A365-D open / A165-D open	on/off battery door; push button	53	128	
		 BTE	A570-D / A370-D / A170-D	on/off battery door; push button; VC toggle	62	132	
		 BTE open	A570-D open / A370-D open / A170-D open	on/off battery door; push button; VC toggle	53	128	

Replacement overview

previous hearing aid family	new hearing aid family	new hearing aid family overview (models available in each Avio series)		controls (standard/option)	gain [dB]	max out-put [dB]	fitting range
STARIS			SA10	on/off battery door	45	121	
		CIC 10; CIC 10 HPG	SA10 HPG	on/off battery door	55	126	
			SA30 / SA30 HPG	on/off battery door; push button; VC wheel	55 / 60	123 / 130	
		ITC 30; ITC 30 HPG	SA30-D / SA30-D HPG	on/off battery door; push button; VC wheel	55 / 60	123 / 130	
D			SA40 / SA40 HPG	on/off battery door; push button; VC wheel	62 / 70	131 / 135	
		ITE HS 40; ITE HF 40 HPG	SA40-D / SA40-D HPG	on/off battery door; push button; VC wheel	62 / 70	131 / 135	
ANALOGUE			SA50 / SA50 HPG	on/off battery door; VC wheel; push button;	62 / 70	131 / 135	
		ITE FS 50; ITE FS 50 HPG	SA50-D / SA50-D HPG	on/off battery door; push button; VC wheel	62 / 70	131 / 135	
			SA70	on/off battery door; push button; VC wheel	71	134	
		BTE	SA70-D	on/off battery door; push button; VC wheel	71	134	
SPRINTER PRO			SA70 open	on/off battery door; push button; VC wheel	62	133	
		BTE open	SA70-D open	on/off battery door; push button; VC wheel	62	133	
		SA80	on/off battery door; push button; VC wheel	72	136		
	P BTE	SA80-D	on/off battery door; push button; VC wheel	72	136		
STARIS/EVO Mega Power			EL90-VI Linear	on/off battery door; push button; VC wheel	86	145	
		SP BTE	EL90-VI WDRC				

Replacement overview

previous hearing aid family	new hearing aid family	new hearing aid family overview (models available in each Avio series)		controls (standard/option)	gain [dB]	max out-put [dB]	fitting range		
ANALOGUE	STAGE	 CIC; CIC HPG	ISG13-P / ISG 113-P	on/off battery door	52	125			
			ISG23-P / ISG 123-P	on/off battery door; push button	52	125			
		 ITC; ITC HPG	ISG33 / ISG 233	on/off battery door; push button; VC wheel	55	125			
			ISG33-P / ISG 233-P	on/off battery door; push button; VC wheel	60	128			
		 ITE; ITE HPG	ISG43-P	on/off battery door; push button; VC wheel	68	136			
			ISG343-P	on/off battery door; push button; VC wheel	68	136			
		SPRINTER PRO	STAGE	 BTE	ISG273-V	on/off battery door; push button; VC wheel	64	131	
					ISG373-V	on/off battery door; push button; VC wheel	64	131	
				 P BTE	ISG283-V	on/off battery door; push button; VC wheel	72	139	
					ISG383-V	on/off battery door; push button; VC wheel	72	139	
		 SP BTE	ISG393-VI	on/off battery door; push button; VC wheel	86	145			
X	INTERTON Bodyworn		IBW2-V	on/off VC wheel; push button; bone conductor oscillator	71	133			
			IBW3-PVI	on/off VC wheel; push button; DAI; LED indicator; bone conductor oscillator	90	150	